



Life's too short to slow down

Our extensive range of scooters and accessories are designed to keep you on the move, for longer.

Q RAA Mobility



Choosing your new set of wheels

When it comes to choosing the right mobility scooter, it's about more than just getting from A to B. It's about comfort, speed and style.

Whether you're considering a heavy duty, mid-sized or portable model, we offer expert advice and free in-home demonstrations to help you make the right choice.

By test driving a few Shoprider models, you'll get a feel for how they operate and which one best suits your needs.

Visit the RAA mobility scooter super site at Mile End, and select metro and country RAA Shops. Did you know we also offer in-home demonstrations? Let us come to you!



Shoprider Rocky 8

Equipped with sleek low-profile tyres for a smooth ride, this scooter features bright LED lights for enhanced visibility in any condition. Its digital display keeps you informed with real-time speed, battery life, and other essential information. With an impressive range, this scooter ensures you can travel further with confidence and comfort, making it the perfect companion for all your mobility needs.



Specifications

Length: 1,590mm
Width: 730mm
Height: 1,285mm
Turning circle: 1,860mm
Wheel size: 330mm

Weight: 147kg
Charger: 8Amp
Motor size: 2Hp
Battery size: 2 x 75Amp
Maximum range: 54km

Maximum user weight flat ground: 225kg
Colour options:
Red, blue, dark grey

Free in-home mobility scooter demonstration

Experience peace of mind as you test the scooter in the comfort of your own surroundings. We have a wide range of Shoprider mobility scooters to suit all needs and budgets.

Available in selected areas.[^]

Call 1800 251 254 to speak with us today



[^]In-home demonstrations available only within the Adelaide metropolitan areas and select regional areas. To find out more about the selected locations or to arrange a free home demonstrations in Adelaide metropolitan area and selected regional areas, call 1800 251 254.

Shoprider Rocky 4

Featuring a fresh, modern design and an ultra-comfortable seat for long rides. Stay connected on the go with a built-in USB port, perfect for charging your devices while you travel. Built for reliability, this scooter delivers consistent performance, ensuring you can depend on it for all your mobility needs.



Specifications

Length: 1,440mm
Width: 700mm
Height: 1,300mm
Turning circle: 1,800mm
Wheel size: 330mm

Weight: 146kg
Charger: 8Amp
Motor size: 2Hp
Battery size: 2 x 75Amp
Maximum range: 54km

Maximum user weight flat ground: 225kg
Colour options:
Red, blue

Shoprider All Rounder

Our mid-size mobility scooter offers the perfect balance of style and functionality, fitting most people comfortably while maintaining a sleek, modern design. Built for dependability, it delivers exceptional reliability, ensuring smooth performance wherever you go.



Specifications

Length: 1,300mm

Width: 640mm

Height: 1,200mm

Turning circle: 1,440mm

Wheel size: 330mm

Weight: 99kg

Charger: 8Amp

Motor size: 1.5Hp

Battery size: 50Amp

Maximum range: 40km

Maximum user weight flat

ground: 159kg

Colour options:

Red, blue, dark grey

Shoprider 889SL

This mobility scooter offers ultimate comfort with a double pillow seat and added privacy with a modesty panel. Powered by a strong motor, it delivers smooth, reliable performance on any terrain – combining comfort, style, and power for your travels.



Specifications

Length: 1,300mm

Width: 640mm

Height: 1,200mm

Turning circle: 1,440mm

Wheel size: 330mm

Weight: 99kg

Charger: 5Amp

Motor size: 1.2Hp

Battery size: 50Amp

Maximum range: 40km

Maximum user weight flat

ground: 136kg

Colour options:

Champagne

Shoprider Explorer

This mid-size mobility scooter is engineered for both reliability and comfort, making it an ideal choice for daily use. All controls are conveniently located on top of the tiller, allowing for easy access and smooth operation while riding. The scooter also features a contoured seat, designed to provide excellent support and comfort during longer rides.



Specifications

Length: 1,280mm

Width: 600mm

Height: 1,070mm

Turning circle: 1,500mm

Wheel size: 260mm

Weight: 99kg

Charger: 5Amp

Motor size: .75Hp

Battery size: 50Amp

Maximum range: 40km

Maximum user weight flat

ground: 136kg

Colour options:

Silver, blue

Shoprider GK10 Crossover

This reliable mobility scooter is designed for easy transport, quickly pulling apart for hassle-free storage and travel. It features a comfortable captain's seat, offering superior support for longer rides, making it the perfect companion for both everyday use and adventures on the go.



Specifications

Length: 1,190mm

Width: 585mm

Height: 960mm

Turning circle: 1,300mm

Wheel size: 260mm

Weight: 78kg

Charger: 5Amp

Motor size: 1.2Hp

Battery size: 36Amp

Maximum range: 27km

Maximum user weight flat

ground: 91kg

Colour options:

Red, blue

Shoprider GK9 Little Ripper

This lightweight mobility scooter is easy to transport, pulling apart effortlessly for convenient storage and travel. Despite its light design, it remains comfortable and reliable, offering a smooth ride and dependable performance wherever you go. Perfect for those seeking mobility with ease and convenience.



Specifications

Length: 1,050mm

Width: 560mm

Height: 880mm

Turning circle: 1,150mm

Wheel size: 230mm solid

Weight: 58kg

Charger: 5Amp

Motor size: 0.8Hp

Battery size: 21Amp

Maximum range: 20km

Maximum user weight flat

ground: 136kg

Colour options:

Red, light blue

Ikon eTravel Scooter

This ultra-lightweight mobility scooter weighs just 17kg without the battery, and is perfect for air travel – folding up easily in one piece for hassle-free transport. Powered by a long-lasting lithium battery, it offers impressive range while maintaining exceptional reliability. Designed for convenience and performance, it's the ideal companion for travellers on the move.



Specifications

Length: 910mm

Width: 420mm

Height: 870mm

Turning circle: 1,420mm

Wheel size: 165/175mm solid

Weight: 19kg

Charger: 2Amp
Lithium-ion supplied

Motor size: 180w

Battery type:
Lithium 9.9Ah

Maximum range: 17km

Maximum user weight flat
ground: 115kg

Colour options: Red, green

Other scooters

Please speak to a consultant for pricing and availability



RainRider

The RainRider has been designed and engineered to maximise the users' comfort and safety whilst still being practical. Reliable, stylish and built for all weather conditions.



Rocky 12

Embrace the independence you deserve with the unparalleled strength and reliability of Shoprider's Rocky 12 scooter. This powerhouse scooter has a safe working load of 350kgs on flat, level ground.



Viking

Combine style and reliability with our 3-wheel mobility scooter, designed to look like a sleek motorcycle. With a powerful motor, stable design, and easy handling, it offers the perfect blend of performance and practicality. Ride confidently and turn heads wherever you go!



Rocky Golf

Don't let mobility challenges hinder your love for golf. Step onto the course with confidence and style, knowing that you have the perfect companion to accompany you on your golfing adventures.

Which model is right for you?

Model	IKON ETravel Folding Scooter	GK9 Littl Ripper	GK10 Crossover	888SE Explorer	Deluxe 889SL	All Rounder	Rocky 4	Rocky 8
Maximum speed (km p/h)	6	6	9	10	10	10	10	10
Portability	✓	✓	✓	✗	✗	✗	✗	✗
Nudge bar	Small	Small	No	✓	X	✓	✓	✓
Seat type	Lower backrest only	Mid back seat	Mid back captain seat	High back contured	Double pillow comfort	High back contured	Large high back contoured	Large high back contoured
LED display	✗	✗	✗	✗	✗	✗	✗	✓
USB charging	✗	✗	✗	✓	✗	✓	✓	✓
Spedometer/odometer	✗	✗	✗	✗	✗	✗	✗	✓
Engine size	180w	0.8HP	1.2HP	0.75HP	1.2HP	1.5HP	2.0HP	2.0HP
Wheel size	165/175 solid	230 solid	260	260	330	330	330	330
Battery type GEL (2)	Lithium 9.9Ah	21amp	36amp	50amp	50amp	50amp	75amp	75amp
Charger	2Amp Lithium-ion supplied	5Amp supplied	5Amp supplied	8Amp supplied	5Amp supplied	8Amp supplied	8Amp supplied	8Amp supplied
Maximum travel range (km)	Up to 17	Up to 20	Up to 27	Up to 32	Up to 40	Up to 40	Up to 40	Up to 54
Maximum user height (cm)	170	160	170	180	190	190	200+	200+
Maximum user weight flat ground (kg)	115	136	136	136	136	159	225	225
Maximum user weight 9 degree gradient (kg)	NA	91	91	91	91	125	159	125
Turning circle (cm)	142	115	130	150	144	144	180	186
Scooter weight (kg)	19	57.6	78	99	99	99	146	147
Overall dimensions	L:910mm x H:870mm x W:420mm	L:1050mm x H:880mm x W:560mm	L:1190mm x H:960mm x W:585mm	L:1280mm x H:1070mm x W:600mm	L:1300mm x H:1200mm x W:640mm	L:1300mm x H:1200mm x W:640mm	L:1430mm x H:1300mm x W:680mm	L:1590mm x H:1285mm x W:730mm

Stay mobile with our range of walking aids

Enduro Powerlite Electric Wheelchair

This carbon fibre powerchair, weighs just 13.2kg and is built for both strength and convenience, supporting users up to 140kg. It features a digital display for easy control and monitoring, and folds up compactly for effortless transport. Reliable and versatile, it can also be operated by an attendant using a Bluetooth remote control, offering flexibility and peace of mind wherever you go.



Specifications

Length: 1,045mm

Width: 555mm

Height: 890mm

Turning circle: 1,080mm

Charger: 2Amp
Lithium-Ion supplied
Motor size: 2 x 250w
Battery type: Lithium
10Ah

Maximum range: 18km
Maximum speed: 6km p/h
Maximum user weight flat ground : 140kg

Redgum Enduro Superlite Walker

Specifications

Carbon Fibre Frame
with handbrake

Wheel size: 200mm

Seat height: 580mm

Overall width: 610mm

Product weight: 5.9kg

Maximum load: 130kg



Redgum Opal Transit Chair

Specifications

Seat height: 580mm

Distance between armrests: 460mm

Product weight: 9.4kg

Maximum load: 136kg

Colour options: Blue, yellow



Redgum Adjustable 8inch Walker

Specifications

Seat height: 530–590mm

Wheel size: 200mm

Overall width: 580mm

Product weight: 6.7kg

Maximum load: 120kg

Colour options: Blue, pewter



Redgum Tri Walker

Specifications

Handle height: 830–890mm

Wheel size: 190mm

Overall width: 675mm

Product weight: 4.8kg

Maximum load: 115kg



Redgum Primo 8inch Walker

Specifications

Seat height: 505–555mm

Wheel size: 200mm

Overall width: 610mm

Product weight: 8.3kg

Maximum load: 136kg

Colour options: Light blue, grey



Ride in safety and style

Shoprider accessories to make your life that little bit easier

Deluxe sun canopy

Enjoy the great outdoors while keeping the sun off your head with this handy canopy.

Note: can only be mounted with the bag carrier.



Bag carrier

The perfect place to store your shopping, be it groceries or a new pair of shoes. The carrier also includes a walking stick holder.



Wheelie walker carrier

Suitable to carry most front folding walkers. Easy to remove and attach.



Safety flag

Be safe and be seen with this safety flag.

Note: can only be mounted with the bag carrier.



Storage covers

Protect your scooter from rain, sun and dust.



Breakdown buddy

Handy 4-in-1 air compressor, charger, light and hazard light.



Premium chargers

Dust and moisture resistant, app connectivity with 5-year warranty.

Help when you need it most

Insurance

Get cover with our Mobility Scooter Insurance so if the unexpected happens, you can rest easy. To get a quote, give us a call on 8202 4567, visit an RAA Shop or search RAA Motorcycle Insurance online.

Issued by RAA on behalf of Allianz South Australia Insurance Limited AFSL 232525. T&Cs, limits and exclusions apply. Consider the Product Disclosure Statements, available with the Target Market Determinations at raa.com.au.

Road Service

If you break down and you hold any level of RAA Road Service, call us and we'll try and get you moving.

If we're not able to get you going again and the breakdown happens away from home, we'll organise a taxi to take you and your scooter to a place of safety.

Terms, conditions and limits apply, as set out in your Road Service Guide.

If your mobility scooter needs servicing or repairs, we're here to help.

Like cars, mobility scooters need to be serviced regularly so they run safely and smoothly.

We offer:

- 6-month minor servicing
- 12-month major servicing
- Replacements only using genuine parts
- Warranty support
- Mobile service (we come to you) or visit our Mile End Service Centre

Plus, members save 15% on servicing and 10% on parts.

**Call us on 8202 4477
to make a service booking**



*T&Cs apply and are available from RAA. Servicing and repairs provided for Shoprider mobility scooters only. Mobile pick-up/drop-off available in metropolitan Adelaide only and mobile servicing subject to the suitability and safety of the location as determined by RAA, additional fees may apply. Warranty support is for Shoprider mobility scooters only.

Let us come to you

Book your free in-home mobility scooter demonstration*



**Call 1800 251 254 to
speak with us today**



Better to be a member

Q RAA Mobility

*In-home demonstrations available only within the Adelaide metropolitan areas and select regional areas.

Notes



Keen to take one for a spin?

Visit the RAA mobility scooter super site at Mile End or select metro and country RAA Shops. For in-home demonstrations and expert advice, call us on **1800 251 254**.